




















Experiment 29

Laser Beam Analysis



1. Introduction	6
2. Experimental Set-up	6
2.1. Set-up with the Diodelaser Module	6
2.1.1. The beam analyser head	6
2.1.2. Diodelaser Module	6
2.1.3. Diodelaser Controller LDC-01	7
2.1.4. Collimating Optics	8
2.1.5. The Beam shaping Optics	8
2.2. Set-up with HeNe-Laser	8
2.2.1. HeNe-Laser and Power Supply	8
2.2.2. Eightfold Beam Expander	8
2.3. Examples of measurements	9
2.3.1. Diodelaser	9
2.3.2. HeNe-Laser	9
3. BeamMaster Introduction	10
3.1. Theory of Operation	10
3.2. Overview	11
3.3. BeamMaster and Beam Diagnostics	11
3.4. Quick Start	11
3.4.1. Quick Start: Complete System	11
3.4.2. Quick Start: Components	12
4. Setup	12
4.1. Hardware Installation	12
4.2. Care & Handling	12
4.2.1. Detector Head	12
4.2.2. Interface Board	12
4.2.3. Filters	12
4.2.4. Laser Safety	12
4.3. Hardware Requirements	12
4.3.1. Installing the Interface Board	13
4.3.2. Configuring Hardware with Windows 95	14
4.3.3. Configuring Hardware with Windows 3.1	14
4.3.4. Connecting the Detector	15
4.3.5. Selecting a Filter	15
5. Software Installation	15
5.1.1. In Windows 3.1:	15
5.1.2. In Windows 95:	15
5.1.3. In DOS:	15
5.2. Running the Software	15
5.3. Configuring Software	16
6. Operation	16
6.1.1. Title Bar	16
6.1.2. Menu Bar	16
6.1.3. Tool Bar	16
6.1.4. Window Area	16
6.1.5. Status Bar	17
6.2. Keyboard Operation	17
6.3. Tool Bar	17
6.4. Window Area	17

6.5.	Status Bar	17
6.6.	Methods of Data Entry	18
6.6.1.	Selecting an Option box	18
6.6.2.	Radio Button	18
6.6.3.	Using the Value Entry Box	18
6.6.4.	Colors	18
7.	Menus and Commands	18
7.1.	File Menu	19
7.2.	 Load Configuration File	19
7.3.	 Save Configuration File	19
7.4.	Configuration File Contents	19
7.5.	 View File	19
7.6.	 Delete File	20
7.7.	Save Snapshot file	20
7.8.	Importing Files into Excel	20
7.9.	Create Attenuator File	20
7.10.	 Print File	21
7.11.	Print Setup	21
7.12.	Exit	21
8.	System Menu	21
8.1.	Low Resolution	22
8.2.	High Resolution	22
8.3.	System Setup	23
8.3.1.	Wavelength	23
8.3.2.	System Clip Level	23
8.3.3.	Average	23
8.3.4.	Run Mode	23
8.3.5.	Beep	23
8.4.	Hardware Setup	24
8.4.1.	Plug-in Address	24
8.4.2.	Interrupt	24
8.5.	Start/Stop Link	25
8.6.	Link Setup	25
8.6.1.	Port	25
8.6.2.	Frequency/Baud Rate	25
8.6.3.	Parity Bit	25
8.6.4.	Data Bit	25
8.6.5.	Stop Bit	26
8.6.6.	Transmit	26
8.6.7.	Making a Null Modem Cable	26
9.	View Menu	27
9.1.	 Power	27
9.1.1.	Attenuator	27
9.1.2.	Offset	27
9.1.3.	Power Setup	29

9.1.4.	Units	29
9.1.5.	Power Range	29
9.1.6.	Alarm Limits	29
9.1.7.	Sound On	30
9.1.8.	Max	30
9.1.9.	Min	30
9.1.10.	Large Display	30
9.2.	 Plot Power	30
9.2.1.	Reset	30
9.2.2.	Save Plot	30
9.2.3.	Plot Setup	31
9.2.4.	Time Scale	31
9.2.5.	Units Plot	31
9.2.6.	Plot Range	31
9.3.	 V & W Profile	31
9.3.1.	Gaussian	32
9.3.2.	Profile Zoom	32
9.3.3.	Profile Setup	32
9.3.4.	Supplementary Clip Levels	33
9.3.5.	Profile Grid	33
9.3.6.	Autoscale	33
9.4.	 Position	33
9.4.1.	Position Setup	34
9.4.2.	Position Zoom	34
9.4.3.	Position Grid	34
9.5.	 Chart Width/Position	34
9.5.1.	Save Chart	34
9.5.2.	Chart Setup	35
9.5.3.	Chart Range	35
9.6.	 Projection	35
9.6.1.	Projection Setup	36
9.6.2.	Dimension	36
9.6.3.	Projection Grid	36
9.6.4.	Axes	36
9.6.5.	Graphic Resolution	36
9.6.6.	Rotation Angle	36
9.7.	Status Bar Command	37
9.8.	Tool Bar Command	37
10.	Options Menu	37
10.1.	 Start Data Collection	37
10.2.	 Stop Data Collection	37
10.3.	 Data Collection Setup	37
10.4.	 Start/Stop Test	38
10.4.1.	Test Setup	38
10.5.	 Step	38
10.6.	Save Settings on Exit	38

11.	<i>Windows Menu</i>	39
12.	<i>Help Menu</i>	39
12.1.	Contents	39
12.2.	Search for Help On	39
12.3.	Active Window	39
12.4.	About	39
13.	 <i>Hardware</i>	39
13.1.1.	BM-7 Silicon Detector Head	39
13.2.	Interface Board	40
14.	 <i>Calculations Summary</i>	41
14.1.	Background Noise	41
14.2.	Relative Total Power	41
14.3.	Beam Width	42
14.4.	Beam Centroid Location	42
14.5.	Gaussian Fit	42
14.5.1.	Gaussian Fit Diameter	42
14.5.2.	Gaussian Correlation Coefficient	42
14.6.	Ellipticity	43
14.7.	Speed Issues	43
15.	 <i>Troubleshooting</i>	43
15.1.	System Errors	43
15.2.	Operating Errors	43
15.3.	Analysis Errors	44

1. Introduction

In almost all cases a Laser application requires a Laser beam with low divergence emitted in the so called fundamental Gaussian mode (TEM_{00}). This is not always guaranteed for each laser and especially for high power Laser systems because the emission may be multimode.

Considering the Laser diodes which possess even elliptical beam shapes it is important to have a tool for analysing the structure of the beam shape for different types of Laser in order to be sure to fulfil the required demands for a particular application. Especially for Laser printing, material processing, fibre optics coupling, optical data storage and Laser pumping the spatial intensity distribution of the laser used must be well defined.

The smallest obtainable spot size diameter of a Laser beam depends beside its wavelength on the purity of the spatial mode spectra of the Laser. In case of the Laser printer the desired resolution will not be obtained if the Laser emits additional transverse modes. If more transverse modes are existing the beam shape is no longer circular and depends on the number of higher transverse modes causing a variety of different patterns to occur. This also means that the intensity distribution is different and the application will not fulfil the desired performance.

Within this project a laser beam profiler based on the knife edge method is used to measure the spatial intensity distribution. As probe Lasers a HeNe-Laser and a diodelaser are used.

The HeNe-Laser emits a pure Gaussian TEM_{00} whereas the diodelaser emits an elliptical beam. By means of different optical elements (beam expander and cylindrical lenses) the influence on the intensity distribution will be measured and discussed.

The collection of the data is performed automatically by a PC with special software. For this reason an extra board is built into the PC which controls the beam analyser, collects and stores the data either on floppy or the hard disk of the PC. Based on the recorded data a 2D/3D graphical presentation of the intensity profile is displayed on the monitor or can be printed out on an ink jet colour printer.

2. Experimental Set-up

The complete experimental system is provided with two Laser sources with various optical elements to expand, collimate, shape and focus the emitted beams.

2.1. Set-up with the Diodelaser Module




Fig. 1: Setup with the Diodelaser Module

2.1.1. The beam analyser head




Fig. 2: Beam Analyser head

The Module MH includes the detector head, computer plug-in card, control software, a 10 % and 0.5 % transmission optical filter. The wavelength range is 400 to 1100 nm and the measurable beam size range ($1/e^2$) is 30 μm up to 9 mm. The unit is mounted on a carrier in such a way that the entrance is in line with the optical axis of the rail system.

2.1.2. Diodelaser Module

The Module A contains a Laser diode that emits 30 mW at a wavelength around 820 nm (*Laser Class 3*). The Laser diode is mounted on a Peltier Element, so that the temperature can be varied from 10-30 °C.